

Wine/Beverage DNA Extraction Kit

#WINEBEV100 100 preps

v1.1 June 2023

Sample: 2mL of wine, beer, cider or other beverage. Volume may be reduced for active ferments where samples are viscous

Format: Spin column

Elution volume: 100 µL

Storage: dry at room temperature (15 - 25°C)

Introduction

This kit is designed to purify DNA from bacteria and yeast from a range of beverages, especially those organisms associated with spoilage. The convenient spin column means the process can be completed within 40 minutes, without the use of hazardous chemicals.

Contents of kit

GT Buffer	30 mL
Proteinase K	11 mg x 2 (add 1.1 mL sterile water or elution buffer)
GBT Buffer	40 mL
W1 Buffer	45 mL
Wash Buffer	25 mL (Add 100 mL reagent grade absolute ethanol)
GS columns (bag)	100
2ml collection tubes	200

Protocol

IMPORTANT NOTES BEFORE USE!

• Pulse centrifuge a proteinase K tube and add 1.1 mL elution buffer to the enzyme. Vortex well and leave for 5 minutes before use. Enzyme may be stored at 4 °C for up to a month, otherwise store at -20 °C.

• Add 100 mL absolute ethanol to the Wash Buffer bottle prior to initial use. Mix and cap tightly.

• Additional requirements: absolute ethanol, microcentrifuge tubes (DNase and RNase-free), pipettes. sterile transfer pipettes.

DNA extraction

- Transfer 2 mL of wine/beer/cider other beverage to a 2 mL microcentrifuge tube. *Reduce this volume to 1 mL for active ferments or high yeast samples (e.g. lees).*
- Centrifuge the samples at 12,000 15,000 x g for 5 minutes.
- Using a disposable 2mL transfer pipette, remove the supernatant and discard.
- Add 200 µL of **GT buffer** to the pellet along with 20 µL reconstituted proteinase K (see notes before use above).
- Vortex well and rake along a microcentrifuge rack to break up the pellet.
- Incubate at 60 °C for ≥10 minutes (see notes below)
- Pulse spin to collect the contents at the bottom of the tube and add 200 μ L of **GBT buffer** to the samples.
- Vortex and incubate at 60 °C for ≥10 minutes. During this incubation, add ~110 µL of Elution Buffer to a microcentrifuge tube and incubate at 60 °C until required at the last step.
- Pulse spin to collect the contents at the bottom of the tube and add 200 μL of absolute ethanol to the samples.
- Vortex samples well. If cellular material is present (e.g. lees material) then centrifuge the samples at 12,000 – 15,000 x g for 30 seconds.
- Assemble a GS column into a 2ml collection tube and transfer the supernatant mixture from the previous step into the top of the GS column
- Centrifuge the column assembly (column plus collection tube) at 12,000 15,000 x g for 1 minute.
- Transfer the GS column to a new collection tube and add 400 µL W1 Buffer to the column.

- Centrifuge the column assembly (column plus collection tube) at 12,000 15,000 x g for 30 seconds.
- Discard the contents of the collection tube and reassemble the column into the tube. Add 600 µL **Wash Buffer** (with ethanol added – see notes before use above).
- Centrifuge the column assembly (column plus collection tube) at 12,000 15,000 x g for 30 seconds and discard the contents of the collection tube as previously.

DNA extraction continued

- Reassemble GS column and collection tube then add 200 μL Wash Buffer to GS column. Centrifuge at 12,000 15,000 x g for 1 minute.
- Transfer the GS column to a microcentrifuge tube and slowly add 100 µL pre-heated Elution Buffer to the centre of the matrix in the column.
- Leave column assemblies in rack on the bench for 3-5 minutes. This step allows the column matrix to absorb the heated elution buffer and release the DNA.
- Centrifuge the column/tube assembly at 12,000 15,000 x g for 30 seconds to elute the DNA from the column.
 Place the tubes in the microcentrifuge so that the caps are facing one direction (opposite to the rotation direction – this prevents cap breakage).
- The eluted DNA is now in the microcentrifuge tube and is ready for PCR or qPCR.

DNA may be stored at 4 $^\circ\text{C}$ for up to 5 days or store at -20 $^\circ\text{C}$ for long term retention.

Notes

- 1. The incubation after adding GT Buffer may be extended to 20-30 minutes, if risk of evaporation is minimal.
- 2. A convenient stopping point is after the addition of the GBT buffer or the absolute ethanol to the sample/GT/GBT mixture. Samples may be left at room temperature overnight.
- 3. Rather than empty contents of collection tubes, GS column may be transferred to a new collection tube (additional tubes available from dnature or other suppliers).