

## Instructions for use (v1.5)

#LYSO-25 LysoFAST DNA isolation reagent, 25mL

#LYSO-BOOST-25 LysoFAST BOOST DNA Isolation reagents, 25mL  
*store at room temperature*

**Please read the notes at the end of the protocol sheet, as well as the prewash notes below before starting your method.**

### ***Prewash***

*If supplied with your reagent, some applications (e.g. bacteria and yeast in some beer and wine samples, as well as plant samples) can benefit from a prewash step using 100-200  $\mu$ L of the prewash buffer.*

*The prewash step is carried out at step 1 or after step 2 (in the applications where samples are first centrifuged).*

- 1. Resuspend your sample in the prewash buffer before repeating step 1 (in the case of bacteria/yeast).*

### **Bacteria/Yeast**

1. Centrifuge 1-2 mL of sample (e.g. 15,000 x g for 5 minutes) or slower (e.g. 3,000 x g for 10 minutes)
2. Remove all supernatant with a fine-tip transfer pipette (prewash may be performed here – see above)
3. Add 80  $\mu$ L ready-to-use LysoFAST reagent (if large pellet evident then add 160)
4. Vortex samples well and leave on bench for 3-5 minutes (see notes at end).

Samples are ready to PCR. Use 2-3  $\mu$ L per 20  $\mu$ L PCR/qPCR

## Tissue / Insects

1. Place sample (small or small pieces) in a PCR tube (0.2 mL) or 1.5 mL microcentrifuge tube
2. Add sufficient LysoFAST to cover the sample (30 – 100  $\mu$ L)
3. Leave on bench for 3-5 minutes (see notes at end). Centrifuge 30 seconds to clarify sample

Samples are ready to PCR. Use 2-3  $\mu$ L per 20  $\mu$ L PCR/qPCR

## DNA cards (plant/blood)

1. Place card punch into PCR tube or 1.5 mL tube (depending on size of punch). Prewash may be performed here – see above
2. Add sufficient LysoFAST to cover the sample (30 – 50  $\mu$ L). Ensure punch is submerged under solution (pulse centrifuge if required)
3. Leave on bench for 3-5 minutes (see notes at end).

Samples are ready to PCR. Use 2-3  $\mu$ L per 20  $\mu$ L PCR/qPCR

## Notes

\* some samples will benefit from a short heating step e.g. 95°C for 1-5 minutes.

For LysoFAST BOOST, samples should be centrifuged at 12-15,000 x g for 30 seconds immediately prior to adding to amplification reaction. Avoid pipetting any of the resulting pellet.

If DNA is discoloured, then dilute 1:5 in 10mM Tris, pH 8 or TE buffer, pH 8.

If large pellets remain, centrifuge at 12-15,000 x g for 30 seconds and transfer supernatant to new tube.

If testing RNA, run cDNA synthesis reaction or one-step RT-qPCR reaction immediately after incubation time. LysoFAST is compatible with PerfeCTa qPCR ToughMix and qScript™ XLT One-Step RT qPCR Toughmix reagents, as well as repliqa HiFi ToughMix for conventional PCR applications.

