




• DETECT • DISCOVER •

## Caffeinate yourself with our GC master mix promotion

Here at dnature we love our Solis Bio Dyne reagents. We love them because they give us great results, are room temperature stable and to be honest, we also like the price. Room temp shipping and live in the fridge for months just adds to their usefulness

The GC master mix:

- is ideal for difficult templates or ones with high secondary structures
- is a complete mastermix
- is room temperature stable so it can stay on your bench for a month (or fridge for 6 !)



*Amplify a broad range of targets with varying GC content. And products up to 3kb !*

## Win a 10 cup coffee card from your favourite barista

It's simple:

- Email [orders@dnature.co.nz](mailto:orders@dnature.co.nz), the code is 04-33-0000S by the way, or order online, for a FREE sample of GC master mix ([dnature.co.nz/gcrich](http://dnature.co.nz/gcrich))
- Run the mix alongside your current PCR reagent and email your comparison results to [support@dnature.co.nz](mailto:support@dnature.co.nz)
- At the same time let us know your favourite barista to collect your coffee card

**We're giving 20 cards away so it's a case of first in, first caffeinated.**

Closes May 30<sup>th</sup> or when 20 cards awarded. Prizes may not be exchanged for cash.

